Technical Specifications		
	CE Version	FDA Version 510(K) Number: K182108
Mode	P1/P2	TENS/MCR
Pulse Duration	50- 400Qs	50- 500Qs
Pulse Frequency	2- 150Hz	1- 120Hz
Intensity	P1: 0- 60mA P2:0- 1000QA	TENS:0- 100mA MCR: 0- 1000QA
Time	1- 60min	1- 60min
Battery Life	More than 4 hours	More than 4 hours

# **MS**tim Arth LGT-235

Portable Electro-Stimulation Therapy Device

1

(Longest ()

### Longest

Guangzhou Longest Science & Technology Co., Ltd. Address: #96 Chuangqiang Road, Ningxi Street, Zengcheng District, Guangzhou, China 511399 Tel: 86-20-66353999 Fax: 86-20-66353920 E-mail: export@longest.cn Website: www.gzlongest.com



Health Life Longest Care



## MstimArth LGT-235

MstimArth LGT-235, works with low frequency electro stimulation (MCR and TENS). It is latest technology pain management stimulator for treatment inflammation and pain with knee arthritis.

## **Applications**

Joint cold

Joint sprain





Tired









#### **Features**

⊘ Both Knees Treated Simultaneously








⊘ Wireless Bluetooth Connection

Quick connection, easy and

#### ⊘ Knee Brace with Comfortable Conductive Fabric



#### ⊘ 2 Modes

TENS -- instant pain management at short time MCR -- long-term use to repair damaged tissue