Certification



Technical Specification

Pressure Range: 1.0~5.0 Bar (Step by 0.1 Bar)
Frequency: 1~22Hz
Vibration handpiece: 1~35Hz
Pre-set Protocol: Yes
Circumference: Yes
Dimension: 39.5 (L) x 39 (W) x 19.4 (H) cm
Standard Weight: 11.3kg
Langages: English, Italian, German, Spanish, French

Longest

Guangzhou Longest Science & Technology Co., Ltd. Address: #96 Chuangqiang Road, Ningxi Street, Zengcheng District, Guangzhou, China 511399

Tel: 86-20-66353999 Fax: 86-20-66353920

E-mail: export@longest.cn Website: www.gzlongest.com



PowerWave LGT-2500S PLUS

The ideal body shaping treatment Acoustic waves provide an ideal body shaping treatment option









Two Channels Up to 5.0 Bar Up to 35Hz

• Applications

- ► Anti-cellulite treatment
- Body reshaping
- ► Skin elasticity improvement
- ► Treatment of skin irregularity after liposuction
- Connective tissue tightening
- ► Scar and wrinkle smoothing
- ► Treatment of stretch marks
- Lymphatic system stimulation










• Effects

- ► Improvement in skin elasticity ► Body contouring ► Treatment of pregnancy stretch marks ► Circumference reduction
- ► Cellulite treatment ► Firming tissue ► Smoothing of scars



PowerWave LGT-2500S PLUS

Acoustic Wave Therapy Device

Acoustic wave treatment uses acoustic sound waves together with the body's own ability to repair tissue by stimulating circulation in all layers of the skin, including blood vessels, connective tissue and muscles.

AWT tightens and forms new connective tissue, reducing cellulite and improving skin appearance.

Benefits

- The acoustic wave is able to massage the tissues in depth with considerable beneficial effects.
- Stimulation of blood circulation and lymphatic
- ► Raising the heat tissue and activation of the elimination of fat cells
- ► Activation of metabolism
- Stimulating the production of collagen, improv-ing elasticity and texture of the skin