### Certification



**Promoting member of ISMST** 



CE (MDD)





ISO13485 / ISO9001



**Health Canada** 



**Therapeutic Goods Administration** 



FDA



**NRTL** 



INMETRO

## **Technical Specification**

Pressure Range: 1.0~5.0 Bar (Step by 0.1 Bar)

Frequency: 1~22Hz Vibration handpiece: 1~35Hz Pre-set Protocol: Yes Circumference: Yes

**Dimension:** 39.5 (L) x 39 (W) x 19.4 (H) cm

Standard Weight: 11.3kg

Langages: English, Italian, German, Spanish, French



Guangzhou Longest Science & Technology Co., Ltd.

Address: #96 Chuangqiang Road, Ningxi Street, Zengcheng District, Guangzhou, China 511399

Tel: 86-20-66353999 Fax: 86-20-66353920

**E-mail:** export@longest.cn **Website:** www.gzlongest.com





**Health Life Longest Care** 







# PowerWave LGT-2500S

#### Acoustic Wave Therapy Device

Acoustic wave treatment uses acoustic sound waves together with the body's own ability to repair tissue by stimulating circulation in all layers of the skin, including blood vessels, connective tissue and muscles.

AWT tightens and forms new connective tissue, reducing cellulite and improving skin appearance.

## **Applications**

- ► Anti-cellulite treatment
- ▶ Body reshaping
- ► Skin elasticity improvement
- ► Treatment of skin irregularity after liposuction
- ▶ Connective tissue tightening
- ► Scar and wrinkle smoothing
- ► Treatment of stretch marks
- ▶ Lymphatic system stimulation











## **⊕** Benefits

- ► The acoustic wave is able to massage the tissues in depth with considerable beneficial effects.
- ► Stimulation of blood circulation and lymphatic
- ► Raising the heat tissue and activation of the elimination of fat cells
- ► Activation of metabolism
- ► Stimulating the production of collagen, improving elasticity and texture of the skin



- ► Body contouring
- ► Cellulite treatment
- ► Firming tissue
- ► Smoothing of scars
- ► Improvement in skin elasticity
- ► Treatment of pregnancy stretch marks
- ► Circumference reduction





