
⊕6 Compression Modes




Technical Specification

Chambers	4
Compression mode	6
Time range	1min - 99min
Pressure range	30 - 150mmHg
Rechargeable battery	6500mmAh (Battery life: more than 5 hours)
Adapter	Medical standard (15V 6A)
Input voltage	AC 100 - 240V 50 - 60Hz
Dimension	270(L)X148(W)X128(H)mm
Weight	1.65KG (without battery)
Interface	2.2" OLED screen
Overlap	Yes
Warranty	Control Unit - 2 years, Garments-1 year

Garment choice and sizes - see separate Garment Folder

Certification





CE(MDD)



Health Canada



ISO13485 / MDSAP



Therapeutic Goods Administration



Guangzhou Longest Science & Technology Co., Ltd.

Address: 5&6F, Building B4, No.11, Kaiyuan Avenue, Science City, Guangzhou Hi-tech Industrial Development Zone, 510530 Guangzhou, China Tel: 86-20-66353999 Fax: 86-20-66353920

E-mail: export@longest.cn Website: www.longest.cn / www.gzlongest.com





Health Life Longest Care

Longest



• Features

- Portable design, easy carry
- Pressure: 30-150 mmHg
- High quality overlap design: 4 chambers
- Real time pressure indication, pressure is adjustabled while running
 Warranty: 2 year for control unit, 1 year for garments
- 6 Compression modes, treatment time 1-99 min

- Device Shell: Fire-proof ABS (94V0 standard)
- Garments surface: 420 Denier Nylon, 0.15mm thickness
- Single cuff use is available
- Easy connect, quite operation, user friendly



• Indications

- Relief muscle pains.
- Increase blood circulation to the treated areas.
- Promote lymphatic fluid movement.
- Remove lactic acid, relieve soreness.
- Improve flexibility, stimulating recovery and enhancing performance.
- Prevent thrombus, reduce limbs drops

• Effect Mechanism















Sleeve









® Product Information



Standard Size Units: mm







