Technical Specifications

Chambers	1-3
Compression mode	4
Time range	1min-24hours or continuous
Pressure range	40-150mmHg
Hold time	0-10s(Step by 1s)
Inteval time	10-60s(Step by 1s)
Interface	5" touch screen
Rechargeable battery	Yes(Lithium battery :11.1V 9000mAh)
Adapter	Yes
Input voltage	AC100-240V 50-60Hz
Dimension	250(L)x184(W)x150(H)mm
Weight	2.7kg

Certification



State Food and Drug Adminstration



CE (MDD)





ISO 13485/ISO9001



INMETRO



Guangzhou Longest Science & Technology Co., Ltd.

Address: #96 Chuangqiang Road, Ningxi Street, Zengcheng District, Guangzhou, China 511399

Tel: 86-20-66353999 **Fax:** 86-20-66353920

E-mail: export@longest.cn **Website:** www.gzlongest.com





Health Life Longest Care

Longest



Introduction of deep vein thrombosis (DVT)

Reason









Deep vein thrombosis is one of the most common complications of bed rest, prevention is urgently needed


Features

- 4 Compression Modes
- 5" Touch Screen
- Reusable & Disposable Sleeves
- Gradient Pressure Setting
- Intermittent Pneumatic Compression
- Pressure of Every Chamber Adjustable
- Rechargeable Battery
- Powered by AC or DC



compression modes with reusable & disposable sleeves meet all your need

⊕Compression Modes



Ser Disposable Sleeve

















*Notice: For one - leg treatment, please use single - head hose.

For Reusable Sleeve











^{*}Notice: For two - leg treatment, please use double - head hose.