

Certification



Promoting member of ISMST



SFDA

Specification


Interface	10.4" LCD touch screen
Pressure Range	1.0 – 4.0 bar
Frequency Range	1 – 17 Hz
Pulse	1-15000
Dimension	52.5(L) × 47.6(W) × 110(H) cm
Standard Weight	27.5 kg
Preset Protocols	Yes
User-defined Protocols	Yes
History Database	Yes
Languages	English & French

Longest[®] Guangzhou Longest Science & Technology Co., Ltd.
Address: #96 Chuangqiang Road, Ningxi Street, Zengcheng District,
Guangzhou, China 511399
Tel: 86-20-66353999 Fax: 86-20-66353920
E-mail: export@longest.cn Website: www.gzlongest.com



PowerShocker LGT-2500F

Non-surgical, cost-effective, powerful & easy-operative technology for soft tissue pain relief.



Health Life Longest Care

Longest[®]

LGT-2500F High-performance Shockwave Therapy Machine

Shockwave Therapy is also known as Radial Pressure Pulse Therapy. LGT-2500F is a compressed air-operated ballistic pulse energy generator transferring pulse energy to the transmitter in the handpiece.

The motion and weight of the projectile are accelerated by compressed air producing kinetic energy that is converted into pulse energy when the projectile strikes an unmoved surface (transmitter). This acoustic pulse energy is transformed into impact energy from the contact area into the therapeutic target tissue radiantly.



Applications

Shoulder

- ✔ Calcifying supraspinatus tendinitis
- ✔ Subscapularis tendinitis (Deltoid bursitis)
- ✔ Frozen shoulder
- ✔ Biceps tendinitis

Lower limbs

- ✔ Patellar tendonitis
- ✔ Shin splints
- ✔ Achilles tendonitis
- ✔ Plantar fasciitis
- ✔ Iliotibial band friction syndrome
- ✔ Anserine bursitis
- ✔ Trigger point

Torso

- ✔ Erector spinae pain
- ✔ Soft tissue pain
- ✔ Trigger point

Upper limbs

- ✔ External humeral epicondylitis (Tennis elbow)
- ✔ Internal humeral epicondylitis (Golf elbow)
- ✔ Finger tendinitis
- ✔ Trigger point

Hip

- ✔ The femoral greater trochanterpain syndrome
- ✔ Piriformis syndrome
- ✔ Sciatica
- ✔ Trigger point

Male reproductive organ

- ✔ Erectile dysfunction
- ✔ Peyronie's disease
- ✔ Chronic prostatitis
- ✔ Benign prostatic hyperplasia

Advantages

- ✔ **Results**
Lasting results for patients seeking no-drug, non-surgical solutions.
- ✔ **High ROI**
In-office treatment offers a more profitable solution without the long downtime of complex surgeries.
- ✔ **Simplicity**
Quick, easy-to-perform, delegable procedure.
- ✔ **Clinical Safety**
Non-invasive, no recovery time, drug-free and painless.

Features

Foldable Touch Screen

10.4" LCD color screen; clear and intuitive display

Rotary Knob

Quick setting and sensitive response facilitate operation

Metal Shell

Fashionable design and high quality

Ergonomic Handpiece

Ergonomic grip and elimination of backward; 2 million shocks tested projectile

Integrated Air Compressor

German KNF technology; silent and stable output

Built-in Air Tank

Ensuring the stable, comfortable and more powerful output

360° Lockable Wheel

Silent and flexible movement; lockable design anchors the device when it's placed

Handlebar

Holding design facilitates movement of the device



Profitable Procedure for Practice Owners

Low Acquisition Cost

1 Year Limited Warranty on Hand Pieces

2 Year Warranty on Console

Low Maintenance Cost

In-office Clinical Support

High Patient Satisfaction

Easy Return on Investment

Why Choose to Lease?

Acquiring clinical devices requires a considerable investment to remain competitive. With a simple and fixed monthly payment plan, a lease can enable a business to match scheduled payments with a practice's anticipated revenue stream.

- ✓ Avoid delays and simplify the complex process associated with the standard loan application process
- ✓ Free capital when the needs of a practice demand cash flow is absolutely essential
- ✓ Structured fixed payments over short or long-term basis
- ✓ Leasing does not impact a business credit line
- ✓ Lease and conserve valuable capital

High profitability - Lease Mode

LGT-2500F supports leasing, and this function accomplishes recirculation of the device. Owners can lend it to end users or clinics to make more profit. Each lease period is 30 days. Costumers can renew a period by inputting the PIN code provided by the supplier.



Accessory

- 1 Handpiece (with 15mm Radial Transmitter)
- 2 20mm Radial Transmitter
- 3 Cleaning Brush
- 4 Metal Pipe
- 5 Screwdriver
- 6 Spare Accessory Kit
- 7 Projectile



Optional Transmitters

Transmitters						
Type	6 A (D = 6 mm Acupuncture)*	15 R (D = 15 mm Radial)	15 F (D = 15 mm Focus)*	15 D (D = 15 mm Deep)*	20 R (D = 20 mm Radial)*	36 R (D = 36 mm Radial)
Penetration Depth	0 – 40 mm	0 – 40 mm	0 – 20 mm	0 – 60 mm	0 – 50 mm	0 – 50 mm

Notice: " * " Stands for Optional Accessories.